



The NNF Guide to

Using 3D Printed Rheometer Tools




Computer rendering and photograph of fractal vane modeled after a Bethe lattice with 24 end points. It sits partially submerged in a cup of material (mayonnaise).

Crystal E. Owens, A. John Hart, Gareth H. McKinley



1. Design choice


Choice of vanes: Vanes with 4-24 arms depending on material needs.





4-arm vane

Sparser structures

- Best for stickier and denser materials
- Easy to insert into material



8-arm vane






24-arm vane

Denser structures

- Best for lower-viscosity materials
- Homogeneous stress field

Stress fields:



Most uniform!



Choice of cups

For larger cups and cups with additional features (like optical windows), a second base mounts to the standard Peltier plate.

Alternative base with textured cup filled with ketchup.




2. About the Parts

They are fabricated by stereolithography 3D printing using a methacrylate-based photopolymer with blue dye.

The material is compatible with acetone, isopropyl alcohol, and water. For best results, rinse but do not soak when washing. Avoid rubbing surface vigorously after it has been soaking in pure solvents.

Keep out of sunlight for long-term use.




Rheometer tools right after printing and before post-processing.

Solvent Compatibility


Percent weight gain over 24 hours for a printed and post-cured (60 °C for 5 minutes)

1 x 1 x 1 cm³ cube immersed in respective solvent:

Solvent	24 hr size gain (%)	24 hr weight gain (%)	Solvent	24 hr size gain (%)	24 hr weight gain (%)
Acetic Acid, 5 %	<1	<1	Hydrogen Peroxide (3 %)	<1	<1
Acetone	<1	2	Isooctane	<1	<1
Isopropyl Alcohol	<1	<1	Mineral Oil, light	<1	<1
Bleach, ~5 % NaOCl	<1	<1	Mineral Oil, heavy	<1	<1
Butyl Acetate	<1	<1	Salt Water (3.5 % NaCl)	<1	<1
Diesel	<1	<1	Sodium hydroxide (0.025 %, pH = 10)	<1	<1
Diethyl glycol monomethyl ether	<1	1	Water	<1	<1
Hydraulic Oil	<1	<1	Xylene	<1	<1
Skydrol 5	<1	1.1	Strong Acid (HCl Conc)	<1	<1



3. Getting Started




Textured Cup


Inner Diameter: $2R_c = 28.5$ mm

Depth: 70 mm


Use with TA's
Concentric Cylinder
Peltier Jacket



- ① Insert in jacket and align screws.



- ② Twist in place, aligning both screws with grooves at the same time.



- ③ The screws will fit **tightly** in the grooves, and do not need to go all the way to the end.



The 8-Arm Fractal Vane

Radius: $R_v = 7.6$ mm

Length: $L = 30$ mm




- ① Insert onto spindle, possibly using force
- ② If required, screw on with threaded rod.
- ③ Calibrate and use as normal.

The specific dimensions of the vane affect performance. An M14 ream should give a good fit. If too loose, it will wobble during experiments, so a tighter fit is better. The tightness will decrease mildly with age if it's 3D printed.



4. Sample Loading

Minimizing end effects and boundary layers with vanes




Insertion procedure:

- For thicker samples, insert slowly to ensure no air pockets are trapped inside vane.
- For much thicker materials, insert to depth of 1000 μm and raise slowly back up to 7500 μm to ensure the interstitial spaces in the vane are filled with sample.



5. After use: Cleaning Tips



Removing the vane


For tight-fitting vanes, if a vane is difficult to remove, it can be gently pried off without damaging the instrument.



Cleaning the vane after tests

The vane may be effectively cleaned with small brushes.

The vane is compatible with acetone, isopropyl alcohol, and water. For best results, rinse but do not soak.



Cleaning the spindle

If required, the spindle may be cleaned after use using a cotton swab wetted with isopropyl alcohol.





6. Data Analysis

Converting Torque-Rotation to Stress-Shear Rate

Variables [and units]:

N number of vane arms	R_v vane radius [m]
L vane length [m]	R_c cup inner radius [m]
Ω rotation [rad/s]	$\dot{\gamma}$ shear rate [s^{-1}]
M Torque [N-m]	σ shear stress [Pa]

Precise conversion:

$$\sigma = M / \left[2\pi R_v^2 L \left[\left(1 - \frac{1.113}{N} \right) + \frac{R_v}{4L} \left(2.75 - \frac{3}{\sqrt{N}} \right) \right] \right]$$

This becomes less accurate for

$$\frac{\sigma}{\sigma_{yield}} > \left(\frac{R_c}{R_v} \right)^2 \sim 3.6$$

$$\dot{\gamma} = \frac{2\Omega}{d \ln T / d \ln \Omega}$$

Recommended to fit $T(\Omega)$ to a polynomial or spline so the derivative is smooth.

Approximate fitting (~ 10 -300% error):

$$\sigma \approx 95,000 M$$

Calculated for 8-arm vane

$$\dot{\gamma} \approx 9 \Omega$$

~ 10 -300% error. Estimated from one Carbopol material.

Key References:

- [1] C. Atkinson, J.D. Sherwood. The torque on a rotating n -bladed vane in a Newtonian fluid or linear elastic medium. Proc. R. Soc. Lond. A (1992).
- [2] C. E. Owens, T. Narayanan, A. J. Hart, G. H. McKinley. Additively manufactured vanes with modified geometries for measurements of yield-stress fluids. Society of Rheology (2018).